

Knowledge base sample: report and diagnosis issues

One of five guides rewritten for customers rather than engineers. Every issue carries a triage tag, which is what lets the assistant decide between a self-serve fix and a handover to a person.

Contents

1. Report Not Received
2. Report Is Taking Longer Than Expected
3. I Need My Report Urgently (ECG)
4. I Need My Report Urgently (Echo)
5. My ECG Report Mentions an Incorrect Chest Lead
6. My ECG Report Mentions a Limb Lead Reversal
7. Doctor Has Asked for a Repeat ECG
8. We Need the Patient's Age & Gender
9. Possible Dextrocardia (Heart on the Right Side)
10. I Want to Understand/Question My ECG Report
11. I Want to Understand/Question My Echo Report
12. I Want to Skip/Cancel a Case

1. Report Not Received

We're On It

The Problem You expected an ECG report, but it hasn't reached you yet.

What You'll Notice No report has come through even though you submitted the ECG a while ago, and you're not sure if it's still being processed or if something went wrong.

What Happens Next / What To Do - Tell us the patient's name and approximate time the ECG was taken — we'll check its status for you. - **If it's still processing:** It's in queue with our medical team and will reach you shortly. No action needed on your side. - **If we can't find any record of it on our end:** - Check that your device had mobile data or Wi-Fi turned on at the time of the recording — a weak or dropped connection is the most common reason an ECG doesn't reach us. - Please resubmit the ECG once you're connected to a stable network. - If you're using a MAC/Communicator device, also check that the device shows good signal strength, is properly connected to the machine, and doesn't have any older unsynced ECGs stuck in its memory. Moving to an area with better network coverage or doing a soft reset on the communicator often helps. - If the ECG still doesn't go through after this, please write to us at support@tricog.com and our technical team will look into it.

2. Report Is Taking Longer Than Expected

We're On It

The Problem Your report is taking longer to arrive than it usually does.

What You'll Notice It's been noticeably longer than your usual turnaround time, and you haven't received an update.

What Happens Next / What To Do - Reach out to us with the patient details and approximate time of the recording. - We'll check exactly where the delay is happening on our end and confirm it's in progress. - If our medical team is experiencing a high volume of cases at that time, we will flag your case for priority reporting. - We'll send you an update once the report has been prioritized, and again once it's ready — you don't need to keep checking in, we'll follow up with you.

3. I Need My Report Urgently (ECG)

Talk To Us

The Problem You have a case that needs to be reported quickly — for example, a patient who needs immediate clinical attention.

What You'll Notice You need the ECG report faster than the usual turnaround time because of the urgency of the patient's condition.

What Happens Next / What To Do - Let us know it's urgent when you reach out — mention the patient details clearly. - We'll immediately acknowledge your request and flag the case with our medical team for priority reporting. - You'll receive an email confirming that your request has been forwarded to the medical team. - We'll keep monitoring the case and inform you as soon as the ECG has been reported by our doctors.

4. I Need My Report Urgently (Echo)

Talk To Us

The Problem You have an Echo study that needs to be reported on priority.

What You'll Notice You need the Echo report faster than usual because of an urgent clinical need.

What Happens Next / What To Do - Let us know it's urgent and share the patient/case details. - We'll immediately acknowledge your request and ask our Echo specialist to prioritize your case. - You'll get an email confirming the case has been forwarded to our medical team for urgent review. - We'll notify you as soon as the Echo report is ready.

5. My ECG Report Mentions an Incorrect Chest Lead

You Can Fix This

The Problem Your ECG report shows a note about an incorrect chest lead placement (for example, V1/V2), but you're confident the leads were placed correctly.

What You'll Notice The report flags a chest lead position issue, even though you placed the electrodes as usual.

What Happens Next / What To Do 1. Double-check the position of the V1/V2 chest leads and reposition them correctly if needed. 2. Check that none of the lead connections are loose. 3. Take the ECG again. 4. If you've checked everything and still disagree with what the report says, tell us — we'll have our reporting doctor personally review it again. - If the doctor agrees a correction is needed, we'll update the report and let you know. - If no change is needed, we'll explain the doctor's reasoning to you.

6. My ECG Report Mentions a Limb Lead Reversal

You Can Fix This

The Problem Your ECG report flags a limb lead reversal, but you believe the leads were connected correctly.

What You'll Notice The report mentions a reversal of the limb leads even though placement looked correct to you.

What Happens Next / What To Do 1. Reposition the limb leads carefully, making sure each one is on the correct limb. 2. Check for loose connections. 3. Try interchanging the leads and re-checking the labels — leads are sometimes mislabeled at the connector. 4. Take the ECG again. 5. If you still disagree with the report after this, let us know — we'll get our reporting doctor to review it. - If a correction is warranted, we'll update the report and inform you. - If not, we'll share the doctor's explanation with you.

7. Doctor Has Asked for a Repeat ECG

You Can Fix This

The Problem Our reporting doctor needs another ECG from the patient before finalizing the diagnosis.

What You'll Notice We'll reach out to you specifically asking for a repeat ECG. This usually happens for one of two reasons:

- **Reason 1 - Signal Quality:** The first ECG had noise or poor quality signals that make it hard to read clearly.
- **Reason 2 - Further Evaluation:** The doctor needs an additional reading to confirm or clarify a finding.

What Happens Next / What To Do 1. If it's a quality issue, make sure the leads are placed correctly and the patient stays still, then take the ECG again. 2. If it's for further evaluation, simply repeat the ECG as guided by us. 3. Submit the new ECG as usual — our medical team will review it and continue the diagnosis from there.

8. We Need the Patient's Age & Gender

We're On It

The Problem Our reporting doctor needs the patient's age and gender to accurately support the diagnosis, and this information wasn't available with the ECG.

What You'll Notice We'll reach out to you specifically asking for this information — it's a routine clinical requirement, not a sign anything went wrong with your ECG.

What Happens Next / What To Do 1. Simply share the patient's correct age and gender with us when asked. 2. We'll pass this on to our medical team right away so they can complete the diagnosis.

9. Possible Dextrocardia (Heart on the Right Side)

We're On It

The Problem Some patients have a rare condition called dextrocardia, where the heart is positioned on the right side of the chest instead of the left. This can make ECGs look unusual across multiple readings even when everything was done correctly.

What You'll Notice If several ECGs for the same patient repeatedly show what looks like limb and chest lead reversal, we may reach out to check for this condition.

What Happens Next / What To Do 1. When we contact you, let us know if the patient has a known history of dextrocardia or any similar heart-position variation. 2. We'll share this information with our medical team so they can read the ECG correctly, accounting for the condition. 3. No further action is needed from you once you've shared this detail.

10. I Want to Understand/Question My ECG Report

Talk To Us

The Problem You have a doubt, concern, or question about something written in your ECG report.

What You'll Notice Something in the report is unclear, unexpected, or you'd simply like more clarity from a doctor.

What Happens Next / What To Do 1. Reach out to us with the ECG ID and a clear description of your concern. 2. We'll gather the details and pass your case on to our medical team for review. 3. If needed, our doctor may request a short call with you to discuss the findings directly — we'll check your availability and schedule it. 4. After the discussion, we'll send you a closure email confirming everything has been clarified. 5. Your case is then marked as resolved.

11. I Want to Understand/Question My Echo Report

Talk To Us

The Problem You have a doubt, concern, or question about something written in your Echo report.

What You'll Notice Something in the report needs more explanation, or you'd like a doctor to walk you through the findings.

What Happens Next / What To Do 1. Reach out to us with the Echo ID/case details and your concern. 2. We'll share this with our medical team for review. 3. Depending on the case, our doctor will either: - Send a written clarification by email, which we'll forward to you for your confirmation, **or** - Request a short call with you to discuss it directly — we'll check your availability and set this up. 4. Once the concern is addressed, we'll send you a closure email. 5. Your case is then marked as resolved.

12. I Want to Skip/Cancel a Case

Talk To Us

The Problem A case cannot move forward for reporting — usually because the uploaded images/study were incomplete and more images can't be provided.

What You'll Notice This usually comes up in one of two ways:

- **Scenario 1:** You uploaded an Echo study, our technician found the images insufficient and asked for more, but the patient is no longer available to provide them (e.g., they've already left).
- **Scenario 2:** You've asked us to generate a report with the images already available, but our technician has determined the images aren't sufficient for an accurate diagnosis.

What Happens Next / What To Do 1. Let us know that additional images cannot be provided. 2. We'll confirm with our medical team that the case cannot be reported with the available information. 3. The case will be closed as **Skipped**, and we'll let you know once this is done. 4. If you're able to arrange the patient for a repeat scan/additional images at any point, just let us know and we'll reopen the case for reporting.

Still Need Help?

If your issue isn't covered above, or you'd like to speak with someone directly, reach out to us at:

support@tricog.com +91 80 4075 3333

We're here to help.